

Product specification sheet



1321-001

SpinChem® RBR S3

Properties

Material:	Stainless steel (SS316L / EN2348)
Diameter:	70 mm
Height:	30 mm
Weight:	397 g (148 g without shaft)
Shaft length:	410 mm
Shaft diameter:	10 mm
Bed volume:	69 mL
Filter porosity:	104 µm

Operating conditions

Rotational direction:	Counterclockwise
Rotational speed:	200 to 800 rpm
NB: Maximum 500 rpm continuously when shaft guide is used!	



SpinChem AB

✉ Tvistevägen 48C, 90736 UMEÅ, Sweden

☎ +46 90192501