

# Product specification sheet

9014-00285



## SpinChem® Essentials kit MiniRBR

1 reaction vessel with a cap made of polypropylene (PP) with 1 MiniRBR

### Tube Vessel

#### Properties

<b>Material Reaction tube</b>	Polypropylene (PP)
<b>Diameter:</b>	0,9 mm
<b>Height:</b>	75mm

### SpinChem MiniRBR

#### Properties

<b>Material:</b>	Glass reinforced polyamide (PA 6, PA 6.6)
<b>Diameter:</b>	16 mm
<b>Height:</b>	20 mm
<b>Weight:</b>	1.8 g
<b>Maximum solids capacity:</b>	1.5 mL
<b>Filter pore size</b>	36 µm
<b>Speed:</b>	300-1200 rpm
<b>Reaction tube working volume:</b>	5-15 mL

#### Operating conditions —

<b>Operating temperature:</b>	-10 to +50 °C
-------------------------------	---------------

