

# Product specification sheet



1321-201

## SpinChem® RBR S3 Plus

### Properties

<b>Material:</b>	Stainless steel (SS316L / EN2348)
<b>Diameter:</b>	70 mm
<b>Height:</b>	63 mm
<b>Weight:</b>	537 g (289 g without shaft)
<b>Shaft length:</b>	410 mm
<b>Shaft diameter:</b>	10 mm
<b>Bed volume:</b>	138 mL
<b>Filter porosity:</b>	104 µm

### Operating conditions

<b>Rotational direction:</b>	Counterclockwise
<b>Rotational speed:</b>	200 to 800 rpm
<b><u>NB:</u> Maximum 500 rpm continuously when shaft guide is used!</b>	



SpinChem AB

✉ Tvistevägen 48C, 90736 UMEÅ, Sweden

☎ +46 90192501