

Product specification sheet



1221-001

SpinChem® RBR S2

Properties

Material:	Stainless steel (SS316L / EN2348)
Diameter:	45 mm
Height:	30 mm
Weight:	317 g (85 g without shaft)
Shaft length:	380 mm
Shaft diameter:	10 mm
Bed volume:	28 mL
Filter porosity:	104 µm

Operating conditions

Rotational direction:	Counterclockwise
Rotational speed:	300 to 1000 rpm
NB: Maximum 500 rpm continuously when shaft guide is used!	

